Prior to use, carefully read the instructions for use.

EN

GC Fuji® II LC CAPSULE

RADIOPAQUE LIGHT CURED REINFORCED GLASS IONOMER RESTORATIVE IN CAPSULES

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- Restoration of Class III, V and limited Class I cavities.
- Restoration of primary teeth.
- Core build up.
- Cases where a radiopaque restoration is required.
- Geriatric applications.
- As a base or liner.

CONTRAINDICATIONS

- Pulp capping.
- Avoid use of this product in patients with known allergies to glass ionomer cement.

DIRECTIONS FOR USE

Powder / Liquid Ratio (g/g)	0.33 / 0.10
Mixing Time (sec.)	10"
Working Time (min., sec.)	3'15"
Light Curing Time (sec.)	20"
Depth of Cure (A2) (mm)	1.8

Test conditions : Temperature (23 +/-1°C)
Relative humidity (50 +/-10%)

1. CAVITY PREPARATION

- Prepare tooth using standard techniques. Extensive mechanical retention is not necessary. For pulp capping, use calcium hydroxide.
- Wash and dry but DO NOT DESICCAT. GC CAVITY CONDITIONER or GC DENTIN CONDITIONER is recommended to remove the smear layer (Fig. 1).

2. CAPSULE ACTIVATION AND MIXING

- Before activation, shake the capsule or tap its side on a hard surface to loosen the powder (Fig. 2).

- To activate the capsule, push the plunger until it is fully depressed with the main body and hold it down for 2 seconds (Fig. 3).

- Note:
- Ensure the plunger is fully pressed to avoid the incorrect mixing ratio of powder and liquid.
 - The capsule should be activated just before mixingand used immediately.

- Immediately set it into a mixer (or an amalgamator) and mix for 10 seconds (+/- 4,000 RPM) (Fig. 4).

- Note:
- When using a rotational mixing device, it is recommended not to use the centrifuge option.

3. RESTORATIVE TECHNIQUE

- Immediately remove the mixed capsule from the mixer and load it into the GC Capsule Applier.
- Make two clicks to prime the capsule then syringe (Fig. 5).

- The working time is 3 minute 15 seconds from start of mixing at 23°C (73.4°F). Higher temperature will shorten working time.

- Note :
- To adjust the direction of the nozzle, hold the applier with the capsule towards you and turn the capsule body.
 - To remove the used capsule, push the applier release button. Twist the capsule and pull upwards.

- Remove surface moisture but DO NOT DESICCAT.

- Extrude cement directly into preparation (Fig.6). Avoid air bubbles.

- Form the contour and place a matrix if required.
- Light-cure for 20 seconds using a suitable visible light curing device (470nm wavelength) (Fig. 7). Place light source as closely as possible to the cement surface.

- Note :
- For cavities deeper than 1.8 mm, use a layering technique.

- FINISHING
Finish under water spray using superfine diamond bur, silicone point and polishing strips (Fig. 8).

- Note :
- Apply a final coat of GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure for 10 sec.) or G-COAT PLUS (light cure for 20 sec.)(Fig. 9).

STORAGE

Recommended for optimal performance, store in a cool and dark place (4-25°C / 39.2-77.0°F).

A1, A2, A3, A3.5, A4, B2, B3, B4, C2, C4, D2 (11 shades based on Vita® shades).

Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

PACKAGES

- Single shade package : 50 capsules.
- Assortment package : 50 capsules (10 pcs each of A4, B2, B3, C2; 5 pcs each of B4, D2).
*Average contents per capsule : 0.33g powder and 0.10g (0.085mL) liquid.
*Minimum net volume of mixed cement per capsule : 0.10mL.

Option : GC Capsule Applier (1 piece).

CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

Undesired effects- Reporting:

If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilance@gc.dental

In this way you will contribute to improve the safety of this product.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.gceurope.com> They can also be obtained from your supplier.

Last revised : 07/2019

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.gceurope.com> They can also be obtained from your supplier.

2797

GC

MANUFACTURED by
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

DISTRIBUTED by
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-KU, Tokyo 174-8585, Japan

EU : GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33,
B-3001 Leuven, Belgium TEL +32 16 74 10 00

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