

For use, carefully read  
the instructions for use.

EN

**G-aenial®**  
LIGHT-CURED RADIOPAQUE COMPOSITE RESTORATIVE

For use only by a dental professional in the recommended  
indications.

**RECOMMENDED INDICATIONS**

A) General:

1. Direct restorative for Class III, IV, V cavities.

2. Restorative use for wedge-shaped defects and root surface  
cavities.

3. Diagnostic for veneers and discoloration.

B) General POSTORSH:

1. Direct restorative for Class I and II cavities.

**CONTRAINDICATIONS**

1. Null pulp.

2. G-aenial® product in patients with known allergies to  
methacrylate monomer or methacrylate polymer.

**DIRECTIONS FOR USE**

1. Shape Selection:

Clean the tooth with peroxide and water. Shade selection should be done prior to use. Select the appropriate G-aenial shades by referring to the G-aenial shade guide.

2. Cavity Preparation:

Prepare cavity using standard techniques. Dry by gently blowing with air.

Note:

For pulp capping use calcium hydroxide.

For bonding: G-aenial and enamel / dentin, use a light-cure bonding system such as G-aenial® Bond.

Follow manufacturer's instructions.

3. Placement of a cavity:

Insert the G-aenial UniTip into the Unip Tip APPLIKATOR equivalent.

Remove the cap and replace.

Note: Use a thin tip.

4. After placement, use a thin tip Unip APPLIKATOR and Unip tip to remove the adhesive from the mouth. This prevents the Unip APPLIKATOR from getting stuck.

5. Maintain pressure on the Unip tip while removing the Unip APPLIKATOR from the mouth. This prevents the Unip APPLIKATOR from getting stuck.

6. Dispensing from a syringe:

After dispensing, use a thin tip Unip instrument. After dispensing, screw syringe plunger anticlockwise by half a turn to release residual pressure inside the syringe. Replace cap immediately after use.

Note: Initially, material can be applied in a single layer to achieve aesthetic restorations using standard shades. For details, refer to the Clinical Hint A.

7. Maintenance:

After removing the adhesive from the mouth, stand for 5 minutes at room temperature.

8. After dispensing, use a thin tip Unip to remove ambient light. Ambient light can shorten the manipulation time.

**Clinical Hints**

A) General:

1. Anterior cavities:

Restorative technique: Use G-aenial UniTip to manipulate the use of one shade alone will be sufficient. In case where a higher degree of transparency is needed, use two shades.

The Outside shade can be selected. See also Examples of Clinical Applications.

B) General POSTORSH:

1. Dispensing from a syringe:

After dispensing, use a thin tip Unip instrument. After dispensing, screw syringe plunger anticlockwise by half a turn to release residual pressure inside the syringe. Replace cap immediately after use.

Note: Initially, material can be applied in a single layer to achieve aesthetic restorations using standard shades. For details, refer to the Clinical Hint A.

2. Maintenance:

After removing the adhesive from the mouth, stand for 5 minutes at room temperature.

3. After dispensing, use a thin tip Unip to remove ambient light. Ambient light can shorten the manipulation time.

4. In the case of small cavities:

Restore using a one shade technique. For small cavities, use G-aenial UniTip to manipulate the use of one shade alone will be sufficient. In case where a higher degree of transparency is needed, use two shades.

The Outside shade can be selected. See also Examples of Clinical Applications.

5. Contouring without Light Curing:

Contour using standard techniques.

Light Cure G-aenial using a light curing unit (Fig. 1). Keep light guide as close as possible to the surface. Refer to the following section "Irradiation Time and Effective Depth of Cure".

6. Dispensing and Polishing:

After dispensing, use a thin tip Unip to remove ambient light, polishing points and discs.

To obtain a high gloss, polishing pads can be used.

**SHADES**

1. 22 shades for anterior:

Standard shade: Xtra-Bright Bleaching White, WH, A1, A2, A3, A4, B1, B2, B3, C1, C2, C3, C4, C5, C6, C7, C8, C9, C10, C11, C12, C13, A1, A2, A3, B1, B2, B3, C1, C2, C3, C4, C5, C6, C7, C8, C9, C10, C11, C12, C13, C14, C15, C16, C17, C18, C19, C20, C21, C22.

2. 6 shades for posterior:

Standard shade: P-A1, P-A2, P-A3, P-A4, P-E1, P-E2, P-E3, P-E4.

Note:

A) All shades are based on "Vita" shade.

\*Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

**STORAGE**

For optimal performance, store in a cool and dark place.

Light Cure G-aenial using a light curing unit (Fig. 1). Keep light guide as close as possible to the surface. Refer to the following section "Irradiation Time and Effective Depth of Cure".

7. Finishing and Polishing:

After dispensing, use a thin tip Unip to remove ambient light, polishing points and discs.

To obtain a high gloss, polishing pads can be used.

**CLEANING AND DISINFECTING:**

MULTI-USE DELIVERY SYSTEM: To avoid cross-contamination between patients this device requires mild-level disinfection.

After use, use a separate device and label for disinfection.

Discard after use.

DO NOT IMMERSE. Thoroughly clean device to prevent drying and accumulation of contaminants.

Do not use autoclave or sterilization cycle.

DO NOT DISINFECT. Use a separate device and label for disinfection.

DO NOT REUSE.

DO NOT DISINFECTION.

DO NOT DISINFECT.

